

SANDA

Industry Leader in Modern Advanced Titration Systems

Proudly Introduces

Precision Dispenser CE2010

Dual Precision Dispenser with Built-In Magnetic Stirrer



Great Features Include:

- ◆ Great for Analytical, Clinical and Testing Laboratories.
- ◆ Two Precision Buret Pumps with Programmable Dispensing Rates.
- ◆ Pumps Can Operate At The Same Time While Using Different Delivery Rates.
- ◆ Delivery Rates Can Be Changed During Titrant Delivery – Great for pH Titrations.
- ◆ Quick and Easy Titrant Changeover, Superior to Syringe-driven Pumps.
- ◆ Front Controls Allow Manual Purging and Filling.
- ◆ Built-in Magnetic Stirrer with Adjustable Speeds.
- ◆ Durable Stainless Steel Cabinet Top.
- ◆ Driven by Simple RS-232 Commands, Almost Any Program Can Control Dispensing.

The precision and versatility of the Precision Dispenser CE2010 has been proven in analytical work with the SANDA F.A.C.T.S. CE2010 Multi-Titrator, and is also compatible with the SANDA FACTS CE2000.

For more information, contact:
SANDA, 4005 Gypsy Lane, Philadelphia, PA 19129
Phone: 800-999-2993, Fax: 215-849-8102; Email: sales@sanda.com
Web Site: www.sanda.com